

INSTRUMENTATION IN NEUROSURGERY (COAGULATORS, DRILLS, CUSA AND RETRACTORS)

Presented by : Noufal Basheer



Coagulators

History

- Thermal energy for haemostasis, dates back to pharaonics,
- The earliest known surgical records - in papyrus documents from Egypt dated as early as 3000 B.C. called **fire drill**—a device which when turned rapidly produced heat along
- Hot iron
- Bovie and Cushing in 1920s
- Lasers in 1960s



Electrocautery

- variety of electrical waveforms
- a constant waveform, - This produces heat very rapidly, to vaporize or cut tissue.
- an intermittent waveform, produce less heat. Instead of tissue vaporization, a coagulum is produced.

Effect of temperature on cells

>40°C	Reversible cell damage
>49°C	Denaturation (irreversible cell damage)
>70°C	Coagulation (collagens convert to gelatin)
>100°C	Desiccation (cells dry as water vaporises, gelatin becomes 'sticky')
>200°C	Carbonisation , eschar formation (pathological 4 th degree burns)


Monopolar Coagulators

- Electrical energy in the range of 250000 to 2 million Hz
- Heating effect Depends upon the density of current
- Size of electrode should be as small as possible
- Fat, bone and air have low water content and hence high resistance
- Ground electrode must have a large area of contact to ensure low current density
- Healing is slower by 2 days, with wound having less tensile strength and larger scar (Vs scalpel cut)
- Increased susceptibility to wound infection

- The active electrode is in the wound.
- The patient return electrode is attached somewhere else on the patient. The function of the patient return electrode is to remove current from the patient safely.

$$\text{BURN} = \frac{\text{CURRENT} \times \text{TIME}}{\text{AREA}}$$

- The current must flow through the patient to the patient return electrode.



Pad Site Location

Choose: Well vascularized muscle mass


Avoid: Vascular insufficiency
Irregular body contours
Bony prominences

Consider: Incision site/prep area
Patient position
Other equipment on patient


Return Electrode Monitoring, actively monitor the amount of impedance at the patient/pad interface and deactivate system

Modalities


cut	Coagulation	blend
Continuous wave	pulsed	Continuous, with resting period
yellow	blue	blue




Endo cut: fractionated cutting under water




Cutting




Contact Coagulation



Loop Excision



Spray Coagulation



Low Thermal Spread/Charring High
← Voltage →


Variables Impacting Tissue Effect

- **Waveform**
- **Power Setting**
- **Size of Electrode**, The smaller the electrode, the higher the current concentration
- **Time**
- **Manipulation of Electrode**
- **Type of Tissue**
- **Eschar**

Safety measures

- **Start up self check**
- **Return electrode continuity monitor**
- **Contact quality monitoring**
- **Return current feedback monitor**
- **High frequency leakage monitor**
- **Earth leakage monitor**
- **Output error monitoring**
- **Smoke filtration**
- **Activation time limit alarm**
- **Do not activate the generator while the active electrode is touching or in close proximity to another metal object**


- Power output should be sufficient to achieve the desired surgical effect but should not be too high, Power requirements vary according to the desired surgical effect, the active electrode size and type




Bipolar Coagulators

- Greater precision and less damage to tissue
- Less power needed
- Current flows through one blade and out through other
- Only the tissue grasped is included in the electrical circuit
- More predictable and less stimulating muscles and nerves
- More effective for coagulating tissue under a layer of fluid
- Radionics vs malis bipolar instruments
sensing device ,no need of irrigation,chances of inadequate coagulation


Modalities




coagulation



cut



ablation



scissors

SS

- Optimum distance between electrodes
- Continuous irrigation with saline
- Charred tissue should be wiped off with moist clothes ,avoid blade to scrape


Tip diameter	use
1.5,2 mm	Large vessels and scalp bleeders,fascia ,muscles
.7-1mm	Dura and brain surface
.5 mm	Tissue close to blood vessels,nerves,and brainstem

Shaft length		
8cm		Brain surface to depth of 2 cm
9.5 cm		Deeper regions
10 cm		TNTS, posterir third ventricle


	micro	macro
Power range	0.1 – 9.9 watts	1-50 watts
adjustability	0.1 watts	1 watt
	Precise point coagulation	Universal use

- Shafts of different length available
- Self irrigating forceps,pre irrigation and post irrigation function
- Jet irrigation systems in haematomas
- Transistorized coagulator system,equipped with themocontrole system(sugita and tsugane)
- Ohta etal,irrigation on when forceps is close
- PTFE coated forceps

FULLY INSULATED TIP




PARTIALLY INSULATED TIP



Bipolar forceps with jet irrigation system



PTFE coated



Complications

- Formation of coagulum
- Adherence of the blood vessel to the tip of the forceps
- Penetration of aneurysm
- Undesirable regional tissue damage due to grounding of current through the body

General principles

- Current flow should not be started till the desired bleeder is reached
- Current setting should be reduced when changing to fine tipped forceps
- Coagulation should be done in a small pool of water
- When irrigated it should not be flooded
- When using on a vessel forceps should be pulsated
- Current should be set as low as possible
- Should be cleaned immediately after use

Argon-Enhanced Coagulation & Cut

- **Decreased smoke, odor**
- **Noncontact in coagulation mode**
- **Decreased blood loss, rebleeding**
- **Decreased tissue damage**
- **Flexible eschar**

Harmonic scalpel

- leading ultrasonic cutting and coagulation surgical device
- using lower temperatures than those used by electro surgery or lasers
- vessels are coapted (tamponaded) and sealed by a protein coagulum
- Coagulation occurs by means of protein denaturation when the blade couples with protein denatures to form a coagulum that seals small coapted vessels



- control of harmonic Scalpels coagulation rate & cutting speed depends on time & force applied to the tissue by the end effector.
- The Harmonic Scalpel uses ultrasonic technology, & energy that allows both cutting & coagulation at the point of impact.
- As compared to electro surgery
 - 1) fewer instrument exchanges are needed
 - 2) less tissue charring and desiccation occur
 - 3) visibility in the surgical field is improved.


The harmonic scalpel uses ultrasonic energy to both cut tissue and coagulate blood vessels




Laser coagulators

- Types :
 - CO₂ laser,
 - Nd : YAG Laser,
 - Argon Laser,
 - KTP laser
- Principle : Photocoagulation
- Explosive tissue vapourisation
- Coagulation, vapourization, haemostasis, Cutting

DRILLS IN NEUROSURGERY



- Development of early tools  trephination to latest motor powered microdrills
- Records of neurosurgery from 3000 BC shows
1st evidence of trephination
hand operated drill in dentistry- 100 AD
- First powered instrument devised by George f. green,
English dentist in 1869
- Sir Heneage Ogilve 1st air powered drill & osteotome
- Robert m hall forest c barber developed modern high
speed drills

Pneumatic high speed drill system

System comprises of

1) Motors



2) Pneumatic control unit with regulator & various connectors

3) Various attachments & dissecting tools

4) Lubricant/diffuser

Mechanism

- Vane type is the hallmark
- Rotor spindle housed in rotor housing
- Vanes are incorporated on lengthwise slots on the rotor spindle
- Speed ranges from 65000 to 100000 rpm
- Speed more than 25000 – bone melts away easily
 - no tactile sensation



Advantages

- Great precision
- Hands are free for the control
- Time saving
- If used properly it is the safest, for both patients & surgeon

Instructions

- Stable body
- Microscope should be positioned in a comfortable operating position
- All loose materials should be removed from the field
- Hand piece should be of light weight & should be held in pen holding position

- Drilling underwater :
 - 1) It allows the neurosurgeon to visualize prospective structures through bone, which becomes semitransparent when adequately hydrated
 - 2) Underwater drilling protects key neuroanatomical structures from thermal injury
 - 3) Irrigation serves to constantly wash the head of the drill bit
- Visualizing critical structures through bone

- **Drilling parallel to underlying structures**

The movement of the drill bit should proceed along the axis of the underlying structure being exposed. the sigmoid, means a predominantly superior-inferior motion, whereas for the middle fossa dura, the motion is in an anterior-posterior plane.

- Drilled part should be in the form of a saucer rather than in the shape of cup

It provides the neurosurgeon with increased visualization


& working angles ,smaller potential space in which a

pseudomeningocele can develop & decreases the sharp

bony edge that may result in skin tightness and possible

wound breakdown.

- Burr should always rotate away from the critical structures
- Choice of drill bit
 - 1) Cutting burrs work more efficiently when removing large amounts of bone
 - 2) Diamond burrs are used
 - when working close to, or potentially close to, critical neurovascular structures.
 - for hemostasis when used briefly without irrigation at a site of bleeding.
 - 3) The size of the drill to use the biggest one the working space safely allows



S t r i g h t



A n g l e d



F
o
o
t
e
d





1

Snap onto motor



2

Click in tool



3

Rotate to lock

Dissecting tools



Acorn,

Cylinder

Ball

Tapered

Match Head

Bone mill



Applications

- Craniotomy
- Correction of craniosynostosis,
- craniofacial anomalies
- Laminectomy, laminoplasty
- Foraminotomy
- removal of osteophytes, iliac crest grafting etc.
- Excision of odontoid in TOO
- Removal of ACP

Complications

- Direct penetrating injury
- Transmission of heat
- magnetic imaging metal artifacts
- Noise pollution
- Transmission of prion diseases

Electric Drill

- More powerful than pneumatic
- Improved overall system weight and balance - cable lighter, more flexible than pneumatic hose
- Reversible direction
- Cable design prevents incorrect connection and assembly



RETRACTORS


Adequate exposure of the target organ represents a laudable prerequisite of every successful operation.

- Hand held
- Self retaining

Hand held retractors


- **Disadvantages:**

- Slipping from the desired position
- Excessive retraction
- Obscuring vision and light
- Inability to maintain in same position for long time




Dural Retractors












Self retaining retractors


- Mechanical retractor mounts for neurosurgery, in 1930s
- Earliest skull mounted system (Demartel, Malis, Heifetz, Edinburgh, Hamby etc.)
 - mounted on burrhole, craniotomy edge
 - Inadequate bone strength, obscuration of the field
- Soft tissue/muscle mounted and pillar and post devices (House and Urban, Weitlaner)
 - less stable, less flexible


- Skull mounted flexbar devices (Dohn and Carton, Apfelbaum) especially useful in Posterior fossa surgery
 - Leyla retractor, Yasargil
 - adjustment difficulties, extreme length of the flexible arms
 - Table mounted flexbar devices
 - modification by Yasargil and Fox
 - Kanshepolsky, U shaped bar
- * head or retractor movement independent of each other


- Headrest mounted flexbar system
sugita, Greenberg,
- Fukushima and Sano ,4 arms on clamp
secured to mayfield headrest

Self Retaining Retractors










Skull mounted Retractor System



Leyla self retaining retractor

- Yasargil
- Self retaining, no assistance needed
- Uniform holding, no pressure irritations
- Upto 5 flexible arms can be used simultaneously
- No obstruction to operative vision
- No restriction of operating area – critical when using microscope

Ball and socket joint for fixing holding rod to operating table





Fixing rod



FF282R coupler



FF273R or FF274R necessary


YASARGIL

Flexible arm



spatula


S
u
g
i
t
a



NEW JERSEY Retractor system

Advantages:

- unique fixation clamp allows unlimited positioning of the retractor arm along the body of the retractor
- Attaches to virtually all self-retaining retractors
- Two retractor blade supports are available ,allowing the use of both flat and round shaft retractor blades
- Provide improved exposure on Posterior Fossa Craniotomies
- Excellent for nerve root retraction during laminectomy procedures



BUDDE Halo Retractor System




Greenberg Retractor System




50-1507 GREENBERG Primary Bar attaches to skull clamp GREENBERG Retractor Kit (50-1500) is supplied with two primary bars.



50-1516 GREENBERG Floating Secondary Bar Attaches to Primary Bar (50-1506) or Secondary Bar (50-1508) or Short Coupler (50-1522). Can also be mounted onto a flat metal edge of up to 1/2" thickness incorporates a universal joint enabling conical adjustment of angle.




50-1524 GREENBERG Large Instrument Holder is a flexbar arm used to support heavy instruments. The cantilever clamp holds instruments with shafts up to 1" in diameter.



50-1509 GREENBERG Long Retractor Arm attaches to floating secondary bars. GREENBERG Retractor Kit (50-1500) is supplied with four long retractor arms.



50-1522 GREENBERG Short Coupler should be used to mount the GREENBERG Handrest (50-1512) to the floating secondary at level.



50-1525 GREENBERG Universal Retractor System X-large Instrument Holder holds instruments up to 2" in diameter.


50-1513 GREENBERG Paddle Tray attaches to long or shortarm (50-1509 or 50-1511)



C
I
O
W
A
R
D



C
l
o
w
a
r
d



Lumbr retraction handle



Caspar retractors






Crank frame Spinal Retractor System

Designed for posterior thoracolumbar surgeries




The BOOKWALTER Retractor.



Provides simple, multi-directional exposure for a broad range of surgical procedures

Tew Cranial/Spinal Retractor



* Capacity to simultaneously manipulate the spinal cord


Brain Retraction Injury

- The incidence of contusion or infarction from overzealous brain retraction is probably 10% in cranial base procedures and 5% in intracranial aneurysm procedures.
- Brain retraction injury is caused by focal pressure (the retractor blade) on the brain leading to
 - 1) Reduction or cessation of local perfusion

Retraction Injury

- Depends upon
 1. shape
 2. number of the retractors
 3. the pressure
 4. duration of the retraction
- The retraction pressures used are usually in the range of 20 to 40 mm Hg
- use of two small retractor blades may provide exposure equivalent to one large blade with a lower retraction pressure

Retraction



Constant pressure

retraction involves readjusting the retractor blade as necessary to keep the pressure constant ,


this type of retraction is naturally suited to retraction pressure

monitoring

Constant exposure

retraction entails setting the retractor blade once without further adjustment.

The brain is allowed to adjust over time to the fixed retractor blade



Goal	Intervention	Comments	
		Advantage	Disadvantage
Anesthesia	Isoflurane Propofol	Well studied Intravenous administration, reliable evoked potential recordings	Hypotensive Hypotensive, expensive
Metabolic depression	Barbiturates Etomidate Hypothermia (mild) ^b	Well studied Rapid recovery Probably quite effective	Delayed recovery ? effective as barbiturates Mildly cumbersome, some delay in extubation
Improve cerebral blood flow	Hypothermia (deep) Mannitol Induced hypertension	Circulatory arrest Effective Effective	Very cumbersome, coagulopathy ? cerebroprotective value Increased bleeding
Minimize Ca ⁺⁺ injury	Nimodipine Hypothermia (mild) ^b	Probably effective Probably quite effective	Hypotension, no intravenous form at present in U.S. Mildly cumbersome, some delay in extubation
Electrophysiological monitoring	Electroencephalogram Evoked potentials	Standardized monitoring Sensitive to cortical and subcortical ischemia	? sensitive to deep ischemia, limits anesthetic options Only certain regions of brain readily monitored

Goal	Intervention	Comments	
		Advantage	Disadvantage
Cerebral blood flow monitoring	Laser-Doppler	Relatively easy to implement	Local cerebral blood flow only, errors likely
	Thermal diffusion	Relatively easy to implement	Local cerebral blood flow only, errors likely
	Transcranial Doppler	Relatively easy to implement	Large vessels only, difficult to monitor ipsilaterally
Increase exposure	Positioning	Practical	Limited application
	Cerebrospinal fluid drainage	Relatively simple	Injury caused by intracranial shifts
	Cranial base removal; osteotomy	Effective	Technically demanding, time consuming
	Hyperventilation	Effective	Decreased cerebral blood flow, alkalosis
Minimize retraction injury	Brain retraction	No injury if appropriate	Safe limits variable
	Brain resection	Effective	Destructive
	Retraction pressure monitoring: <40 mm Hg	Effective	Mildly cumbersome
	Intermittent retraction: <15 min on, 5 min off	Probably effective	Delays surgery, ? increased reperfusion injury
	Multiple retractors	Possibly effective	More cumbersome than fewer retractors

ULTRASONIC SURGICAL ASPIRATORS



- The original ultrasonic aspirator was developed in 1947 for the removal of dental plaques.
- field of eye surgery in 1967 ,based on the principle of phaco-emulsification.
- First developed in 1976 in the US
- Suction device with a tip that vibrates at ultrasonic speed
- Sonic energy disrupts and fragments

- A console and handpiece
- Console has the ultrasonic generator- 2 types

electrostriction,	magnetostriction
piezoelectric ceramic crystals	change in dimensions of a magnetostrictive transducer
crystals decay	not subject to decay

- Titanium tip vibrates longitudinally at a speed of 23 to 35 khz ,
- amplitude of 100 – 300 microns ,function of setting the vibration level
- small amplitude + disruptive effect restricted to



- Hand piece,
stright vs Angled
short Vs Long
internal vs external coaxial irrigation system
different frequencies
- Irrigation system
to suspend the fragmentaed tissue, to
cool the transducer and to prevent the blockage
of suction system

Mechanism

- Simultaneously fragment, emulsify and aspirate parenchymal tissue rapidly
- Vacuum effect
- Cavitation
- Rupture
- Susceptibility depends upon-
water content
sensitivity to vibration

Fat and brain easily disrupts Vs vessel and nerves

- Tissues with weak intracellular bonds, such as tumors and lipomas, are easy to fragment, whereas tissues with strong intracellular bonds, such as nerves and vessel walls, are difficult to fragment

- Low frequency  high amplitude
Useful in hard and partially calcified tumors
 - High frequency  low amplitude
useful while working near vital structure
- adjustments of the vibration energy, irrigation rates and the suction pressures along with the use of appropriate hand piece optimizes the use

Advantages

- Minimizes,
 mechanical manipulation
 Traction on adjacent tissue
- Avoids thermal injury of cautery
- Clear and less crowded operative field
- Vs laser UA are faster ,good visualization of tumor brain interface. Laser is more precise
- Suitable for HPE as they are not significantly distorted

Complications

- Penetrating injury
- ? Transmission of ultrasonic energy to adjacent vital structures through bone
- Reports of multiple cranial nerve palsies

“Winning is overrated. The only time it is really important is in surgery and war.”

Thank you